FAB Titanium G5

DE - Gebrauchsanweisung FAB Titanium G5	2
EN - Instruction for use FAB Titanium G5	3
FR - Mode d'emploi FAB Titanium G5.....	3
CZ - Návod k použití FAB Titanium G5	4
EL - Οδηγίες χρήσης FAB Titanium G5	5
ES - Instrucciones de uso de FAB Titanium G5.....	6
HR - Upute za uporabu proizvoda FAB Titanium G5.....	7
HU - Használati útmutató a FAB Titanium G5 termékhez.....	8
IT - Istruzioni per l'uso di FAB Titanium G5.....	9
PL - Instrukcja stosowania stopu FAB Titanium G5.....	10
PT - Instruções de utilização do FAB Titanium G5	11
RO - Instrucţiune de utilizare FAB Titanium G5	12
SE - Bruksanvisning för FAB Titanium G5.....	13
DA - Brugsanvisning til FAB Titanium G5	14
ET - Kasutusjuhend FAB Titanium G5.....	15
LT – Naudojimo instrukcija FAB Titanium G5.....	16
LV - Lietošanas pamācība FAB Titanium G5.....	17
NL - Gebruiksaanwijzing FAB Titanium G5.....	18
SK – Návod na použitie zliatiny FAB Titanium G5	19
SL – Navodila za uporabo FAB Titanium G5	20
TR - FAB Titanium G5 Kullanım Talimatı	21

MINDFAB

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EN - Instruction for use FAB Titanium G5

Intended use
FAB Titanium G5 is a medical device for the manufacture of crowns and bridges.

Only for professional user (Dental Technician, Dentist).

The intended patient group provides for persons with partially or non-dentate jaws.

Indication
Crowns and bridges for metal-ceramics up to 16 units (max. 4 units in the anterior region and max. 3 units in the posterior region), superstructures, bars.

Contraindication
In case of known allergic reactions to any of the ingredients.

Milling
FAB Titanium G5 is construed for CNC milling machines. Please follow the instructions and parameters of the respective manufacturer of CAM Software and CNC milling machine.

CAD
The design should be done with appropriate CAD software. Please consider an anatomically reduced framework design for the veneering with ceramic. The wall thickness should not be less than 0.3 mm. Choose a sufficient connector dimension (6-9 mm²). Sharp edges and undercuts should be avoided.

Cutting out the frameworks from the Blank
Remove the milled frameworks with suitable cutting tools and smoothing the supports.

Laser welding
Sufficient and secure connections can only be done by laser welding. The use of Ti-welding wire is an option.

Preparation before ceramic veneering
The frameworks can be elaborated with standard carbide cutters, look for smooth transitions and avoid overlapping material. Please use the same cutter for one alloy to avoid contamination. The minimum thickness of the prepared coping should not be less than 0.3 mm. It's recommended to sandblast the frames with 110 µm of Aluminium oxide with 2-3 bar and clean with steam cleaner. Consider the firing process with max. 800°C, otherwise lattice transformation might occur.

Handling conditions / Safety
Metal dust is harmful to health. Use when grinding and sandblasting dust extraction and respirator with filter FFP3 – DIN EN 149.

Residual risks and side effects
If the instructions are observed during the production processes, incompatibilities with Titan alloys are extremely rare. In case of a proven allergy against an ingredient of this alloy, the alloy must not be used for safety reasons. In exceptional cases, electrochemically induced, local irritations have been reported. When different alloy groups are used, galvanic effects might occur. Please inform your dentist regarding the residual risks and side effects. Any serious incident that involves the product must be reported to the manufacturer and the competent authority in the accorded country. The SSCP is available at <https://ec.europa.eu/tools/eudamed>.

Disinfection of the dental prosthesis before insertion
Workpieces from the dental laboratory must be subjected to immersion or spray disinfection before insertion into the patient's oral cavity and then rinsed under running water.

Single-use
Used discs should be not processed for further use as a medical device.

Disposal Instructions
Please dispose of metal residues and dust in an environmentally friendly manner. Do not allow waste to enter groundwater, water or sewage systems. Contact waste exchanges for recycling. Outer packaging can be disposed of in paper waste.

Storage conditions
Temperature, humidity or light has no effect on the product properties.

Our information and recommendation are based on the state of the art in science and technology and has to be considered correct to the best of our knowledge and experience on this day. The above version shall replace any previous versions.

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Adentatec GmbH
MINDFAB GmbH